



# Super EBA®

**SuperEBA® is an all-purpose, alumina-fortified material for the final cementing of a single crown or inlay. It can also be used as a base, liner, or temporary dressing for pain relief.**

## About Super EBA

SuperEBA's neutral pH and mild formulation provide additional patient comfort during the healing process. Its hydrophilic nature and low film thickness provide unsurpassed convenience, making it an excellent choice for sensitive patients. It adheres to tooth structure in the presence of moisture. Not to be used for implants.



## Tips

Regular set has a slab set time of over 15 minutes, and a shorter in-the-mouth set time of 10-15 minutes. The fast set has a 10-15 minute slab set, and a faster in the mouth set of 4-8 minutes. Film thickness 25 microns.

| PRODUCT                        | QUANTITY | ITEM NUMBER |
|--------------------------------|----------|-------------|
| Standard Kit                   |          | 0921007     |
| Standard Kit (All Fast Set)    |          | 0921007F    |
| Standard Kit (All Regular Set) |          | 0921007R    |
| Powder Fast Set                | 15g      | 0921008     |
| Powder Regular Set             | 15g      | 0921009     |
| Liquid Only                    | 6.5ml    | 0921010     |
| Powder Fast Set                | 2 oz     | 0921012     |
| Powder Regular Set             | 2 oz     | 0921013     |

**keystone**  
industries

[dental.keystoneindustries.com](http://dental.keystoneindustries.com)