

Tecnoflask

Microwavable Flask

Tecnoflask allows Diamond D® and other acrylics to be processed in much quicker time inside of a microwave.

The TF-100 muffle will give your laboratory a cleaner environment and make the process more enjoyable for the production of your work. Included with the kit and available separately is the PL-100 protection plate, which extends the muffle's service life with extra needed protection.



FEATURES AND BENEFITS

- Acrylic material designed for use in microwave
- Cannot be damaged by exertion of pressure
- Polymerization time shortened to four minutes
- Used in the making of total upper & lower dental prostheses
- No isolating material needed between the muffle & the plaster
- Does not require a specific resin

ORDERING INFORMATION

Tecnoflask Kit	#1009198
Tecnoflask Muffle TF-100	#1009200
Tecnoflask Press Board PL-100	#1009199


RECOMMENDED FOR USE WITH



Diamond D Denture Acrylic



Diamond D Ultra Sep



keystone
industries

US: 1.800.333.3131 | EU: +49 77 31 91 21 01 | keystoneindustries.com